Sample Configurations on Exome Sequencing Meghan Towne, MS, CGC; Melissa Holman, MS, CGC; Andrew Giles, MS, CGC; Brooklynn Gasser, MS; Ginger Tsai, MS, CGC; Brian Schoenfeld, MS, CGC Ambry Genetics, Aliso Viejo, CA mtowne@ambrygen.com **TAKE HOME POINTS** BACKGROUND Exome sequencing (ES) requires assessment of hundreds of potentially Ongoing assessment of new evidence can aid in closing the gap for relevant variants diagnostic differences between FSC Trio-based ES that involves complete sequencing of both biological Evidence-driven reclassification increased diagnostic yield for all FSC parents ("parental trio") increases the diagnostic potential of ES Primary drivers were new GDRs and familial co-segregation Real-time inheritance information may clarify variants of uncertain studies significance (VUS) Collection of parental trios can be challenging for a variety of reasons • >25% were due to other evidence sources, underscoring the importance of ongoing, dynamic variant assessment

Me

0

ヒ

Optimizing Diagnostic Potential: Impact of Familial

Aims: Investigate how different familial sample configurations (FSC) impact the clinical utility and reclassification of ES

- When available, parental trios maximize the clinical utility of ES
- 30% of ES cases in this cohort were not a parental trio

METHODS 10,921 ES between 2011 - 2021 Non-parental Parental Trio **Proband-only** Duo trio SC 11% 71% 13% 5% L (n=1158) (n=7717) (n=1451) (n=595) Current Original New evidence supporting ES ES Result (+/-ES Result reclassification reclassification) Analyses for this study

FIGURE 1: DIAGNOSTIC RATES BY FSC

Ambry Genetics[®]

∆ +2%

19%

Original

21%

Current

Proband-only



Negative rates were similar across FSC, suggesting appropriate variant reporting, but a lack of segregation evidence leads to more VUS

FIGURE 2: CHANGES IN DIAGNOSTIC RATES BY FSC

∆ +3%

19%

Current

Duo





